

Logic Functions for Coincidence Logic Unit

User Manual

Version 1.5



Universal Quantum Devices

295 Hagey Boulevard, 1st Floor, West Entrance

Waterloo, ON, Canada N2L 3G1

Phone: +1 519 404 6844 Email: info@uqdevices.com Url: www.uqdevices.com

Empowered by

`dotfast-consulting.at`

Thomas Lehner, Dotfast Consulting

Table of Contents

Table of Contents.....	2
Revision History	2
Timetag Explorer Interface	3
Logic Tab	3
Aprox. Cycle time	3
Coincidence Window.....	3
Coincidence Patterns	3
Delay Tab.....	4
Software Interface	4
Initialization	4
Counter Read Out.....	5
Outputs.....	6
LabView Implementation of logic functions	7
Thresholds	10
Delays	10
Counts	10
Logic Pattern	10
STOP	10
Other settings	10


Revision History

Date	Document Version	FPGA Version	DLL Version	Remark
14.08.10	1.0			Initial revision
19.05.11	1.1		2.14	Put logic functions into separate dll Definition of SetOutputPattern() changed
26.1.12	1.2			Defined output delay
22.3.12	1.3			Add information about Timetag Explorer
09.10.13	1.4	2.17	2.17	Update information about new drivers
19.02.14	1.5	2.17.3	2.17.1	Add description of standard implementation and LabView and Visual Basic code.

Timetag Explorer Interface

Logic Tab

This tab is only available with the UQD devices that include the special “Logic” feature. It allows testing and studying of some of the coincidence features of the logic, like measuring single and coincident detections from arbitrary channels, and adjusting delays between channels. The output on the right side gives the single count rates per cycle time and in thousands per second (kHz) for each Input channel and, below, the count rates for the arbitrary coincidence Patterns.



Aprox. Cycle time

The software tries to read the data with this measurement time. Because of computer processing time, the actual cycle time is different. This time *Cycle / ms* is displayed on the right and changes slightly. Only the actual cycle time is used to calculate detection rates.

Coincidence Window

Here the width of the coincidence window is defined. This is the maximum time from the first detection to the last detection in the coincidence event. When all logically selected detections occur within that time window, then the event is counted. Otherwise the event is not counted.

Coincidence Patterns


Each line of buttons defines one coincidence event. Clicking on the button changes the role of the corresponding input.

+ The input must have an active edge in the coincidence window for the event to be counted.

- The input must not have an active edge in the coincidence window for the event to be counted.

Delay Tab

This tab allows the user to define individual delays for each input channel.



This feature is very useful in the logic mode, where the delays of individual channels can be fine-tuned to maximize coincidence signals.

Software Interface

The interface functions are located in the dll `ttInterface.dll`. There is also support for C++, both for Windows and Linux environments.

A normal implementation of the Logic mode proceeds as follows:

Create a `TTInterface` object, open the `TTInterface` using method `Open()`, then create a `Logic` object and call method `SwitchLogicMode()`. Next use the `TTInterface` object to get the time resolution of the Coincidence Logic Unit with `GetResolution()` and set the count thresholds for each channel with `SetInputThreshold()`. Back in the `Logic` object, set the delays for each channel with `SetDelay(int input, int delay)` and set the coincidence window with `SetWindowWidth(int window)`, both in units of the time resolution. The logical pattern of counts to be recorded should be set up also (see description under `CalcCount`). Now, in a repeating loop, call `ReadLogic()` to capture data and `CalcCount(int pos, int neg)` or `CalcCountPos(int pos)` to return the number of counts for each logical pattern. The loop can be delayed to count for a certain time interval, and use `GetTimeCounter()` to find the actual counting time for each loop, which will be slightly longer than the set delay due to computer processing time. After the counting loop is complete, close the `TTInterface` object with `Close()` to end the session.

Initialization

```
public Logic(ttInterface t)
```

The constructor takes a `ttInterface` object which must be already created, and creates a `Logic` object. The `ttInterface` object can already be opened or not.

```
public void SwitchLogicMode()
```

This method must be called before the other logic functions can be used.

```
public void SetWindowWidth(int window)
```

This method sets width of the coincidence window. The parameter `window` is given in internal units (156.25 ps). The window width can be set as high as $2^{24}-1$.

```
public void SetDelay(int input, int delay)
```

This method sets the virtual delay line for each channel. The parameter `delay` is given in internal units (156.25 ps).

The parameter `input` identifies the channel and ranges from 1 to 16.

Counter Read Out

```
public int[] ReadLogic()
```

This method reads the counter out of the device. The return parameter is not needed (used for debug purposes only).

The data is automatically stored in the `Logic` object for later processing.

This method can be called any time after calling `SwitchLogicMode()`. The time between calls to `ReadLogic()` defines the measurement time interval of the captured data.



The device uses a double buffered memory system. The new measurement starts immediately and is running in the background during the call to `ReadLogic()`. For this reason not a single pulse is omitted during readout.

```
public int CalcCount(int pos, int neg)
```

This method operates on the data read by the last call to `ReadLogic()`. It calculates how often a certain event pattern has occurred in the last measurement interval.

The parameters `pos` and `neg` are bit coded. The rightmost bit corresponds to input 1, the next to input 2, and so on.

The parameter `pos` is given by the integer value of the binary number defined by the bit coded Coincidence Event. It defines the inputs that must have an active edge in the coincidence window for the event to be counted, with range [0..65535].

The parameter `neg` defines the inputs that must have no active edge in the coincidence window for the pattern to be counted. The parameter `neg` is optional, with range [0..65535].

pos Coincidence Pattern	pos	neg	Coincidence Event
0000000000000001	1	0	Singles of input 1
0000000000000010	2	0	Singles of input 2
0000000000000100	4	0	Singles of input 3
0000000000001000	8	0	Singles of input 4
0000000000000011	3	0	Coincidence of input 1 and 2
0000000000000111	7	0	Coincidence of inputs 1, 2, and 3
0000000000000011	3	4	Coincidence of inputs 1 and 2 without interaction of input 3

```
public int CalcCountPos(int pos)
```

This method equals `CalcCount` except that the `neg` parameter is always 0. It has better run time.

```
public int GetTimeCounter()
```

The time counter measures the time between the last two calls to `ReadLogic`. The result is given in multiples of 5 ns.

The time when `ReadLogic()` is called will vary over time because of the limited realtime performance of personal computers. The correct count rates can be obtained when the count values from `CalcCount` are divided by the result from `GetTimeCounter()`.

Outputs

```
public void SetOutputWidth(int width)
```

This function defines the length of the generated output pulses. The parameter is given in 5 ns increments. The maximum value is 255.

```
public void SetOutputPattern(int output, int pos, int neg)
```

This function sets the pattern of the output pulses based on the input coincidence events.

The parameter `output` identifies which output to change, with range [1..3].

The parameter `pos` is the bit coded value of the Coincidence Event of the inputs that must be present, with range [0..65535] as for `CalcCount`.

The parameter `neg` is the bit coded value of the Coincidence Event of the inputs that must *not* be present, with range [0..65535].

Method call	Result
<code>SetOutputPattern(1, 1, 0)</code>	Output 1 fires when input 1 is present
<code>SetOutputPattern(2, 6, 0)</code>	Output 2 fires when inputs 2 and 3 are present within the coincidence window
<code>SetOutputPattern(3, 6, 8)</code>	Output 3 fires when inputs 2 and 3 are present and input 4 is not present at the same time

```
public void SetOutputEventCount(int events)
```

This function is an advanced feature that can be used to fine tune the delay of the outputs. Roughly speaking, `events` is the number of events you expect to occur in one 5 ns time slice under worst-case conditions.

Increasing the value by one will increase the delay of the outputs by 10 ns. With the standard setting of 5 the output delay will be 350 ns.

When the input rate is too high for the given event count, the `OutTooLate` Error Flag will be raised. No output pulse will be generated in this condition. There is never an output pulse generated in the wrong timing. For this reason `OutTooLate` is more a kind of warning, not a hard error condition. It just means that not all pulses are generated.

Hint: When two events occur at the exact same time, the event coming from the input with the smaller input number will be processed first. For this reason input 1 has a slightly better real time performance than input 16.

Minimal implementations of logic functions

These implementations consist of a small but useful script that opens the Coincidence Logic Unit, sets some parameters, and reads once the single and coincidence counts on channels one and two.

Visual Basic

```
Imports TimeTag
```

```
Public Class Form1
```

```
    Dim myTagger As TTInterface
    Dim mylogic As Logic
    Dim Resolution As Double
    Dim Errors As Integer
```

```
    Private Sub Form1_Load(sender As Object, e As EventArgs) Handles MyBase.Load
```

```
        Dim Singles1, Singles2, Coincidences As Integer
        Dim FPGA_v, TimerCounter As Double
```

```
        'open connection to device
        myTagger = New TTInterface
        myTagger.Open()
```

```

'get some basic info and display
Resolution = myTagger.GetResolution()
FPGA_v = myTagger.GetFpgaVersion()
MsgBox("FPGA Version: " & FPGA_v.ToString & vbCrLf & _
      "Timing Resolution: " & Resolution.ToString & " s" & vbCrLf
      & vbCrLf & vbCrLf & _
      vbTab & vbTab & "Start Measurement Now?")

'set input thresholds for channel 1 and 2
myTagger.SetInputThreshold(1, 0.3)
myTagger.SetInputThreshold(2, 0.3)

'initiate logic mode
mylogic = New Logic(myTagger)
mylogic.SwitchLogicMode()

'conicidence window (in units of 156ps)
mylogic.SetWindowWidth(20)

'count events for 1 sec:
'take a first snapshot of counters
mylogic.ReadLogic()

'delay thread for about 1sec
Threading.Thread.Sleep(1000)

'take second snapshot of counters
mylogic.ReadLogic()

'how long was the true measurement time between the two ReadLogic
  comands? Given in multiples of 5ns, therefore we multiply by
  this unit to obtain sec
TimerCounter = mylogic.GetTimeCounter() * 0.000000005

'retrieve observed counts for singles 1,2 and coincidences 12
singles1 = mylogic.CalcCountPos(1)
singles2 = mylogic.CalcCountPos(2)
coincidences = mylogic.CalcCountPos(3)

'show data
MsgBox("Measurement Time: " & vbTab & vbTab &
TimerCounter.ToString("F4") & " s" & vbCrLf & _
      "Singles Channel 1:" & vbTab & vbTab & Singles1.ToString &
      vbCrLf & _
      "Singles Channel 2:" & vbTab & vbTab & Singles2.ToString &
      vbCrLf & _
      "Coincidences Chs 1&2: " & vbTab & Coincidences.ToString)

Close()

End Sub

Private Sub Form1_FormClosing(sender As Object, e As
FormClosingEventArgs) Handles Me.FormClosing
  MsgBox("Thanks for testing the counter. Bye")



  ' close connection to device
  myTagger.Close()

End Sub

End Class


```

LabView



Multi-Channel LabView Implementation of logic functions

This description is based on the LabView VI **logic-ver.vi** (current *ver* is 217) and sub-VI **generate logic pattern2.vi**. Below is the standard front panel.



Thresholds

Thresholds allows the voltage threshold to be set for each channel. Two standard voltages are NIM and TTL.

Delays

Delays allows the delay for each channel in nanoseconds to be set, in order to align coincident counts between channels. The delays can only be set to multiples of the Coincidence Logic Unit's resolution.

Counts

Counts displays the number of events in the last Measurement Time corresponding to a given logic pattern.

Logic Pattern

Logic Pattern is the pattern of coincidences to report in counts. Channel 1 is the leftmost button and channel 6 is the rightmost (expand the control to see more channels). Green means an event on the channel must be present, red means an event on the channel must not be present, and black means to ignore the channel. For example, the first three patterns in the above screenshot are for single counts on channel 1, single counts on channel 2, and coincident counts between channels 1 and 2, respectively.

STOP

Always use the **STOP** button, rather than LabView's Abort Execution so that the connection to the Coincidence Panel Logic Unit is closed correctly.

Other settings

One Measurement Only, when enabled, returns the counts for only one Measurement Time, and then stops execution (otherwise measurements run continuously). This is useful, for example, when changing settings between each measurement.

Timetagger Resolution reports the resolution of the Coincidence Logic Unit in nanoseconds.

Coinc. Window controls the time window in nanoseconds within which counts on different channels are counted as coincident.

Actual Coinc. Window reports the coincidence window actually used in nanoseconds, which must be a multiple of the Coincidence Logic Unit's resolution.

Measurement Time chooses the time in seconds between reading the Coincidence Logic Unit's buffer, meaning Counts are displayed as summed over the Measurement Time.

Actual Measurement Time reports the real time in seconds between reading the Counts, due to finite processing time of the PC.