

Universal Quantum Devices
295 Hagey Boulevard, 1st Floor, West Entrance
Waterloo N2L 3G1
ON Canada
Phone: +1 519 404 6844
Email: info@uqdevices.com
Url: www.uqdevices.com



Our ref.

V1.1 October 2013

Quickstart guide for UQD-Logic-16

How to get started:

1. Turn off the power of the NIM crate, and insert the logic unit.
2. Connect the unit to the computer using the USB cable.
3. Power up the NIM crate. Windows should recognize the new device and configure the drivers.
4. If Windows does not recognize USB device, please unzip "WinUSB Driver" to a folder, chose "manual driver selection" and manually select the folder with "WinUSB Driver". Windows will then configure device.
5. Install the software "TimeTagExplorer" to test and operate the unit directly.
6. Make sure the latest version of "Microsoft DotNet Framework" is installed on the Windows system.

